



Now available in
an XR range


PROCLEAR® MULTIFOCAL. ALL-DAY COMFORT FOR EVEN MORE OF YOUR PRESBYOPIC PATIENTS.

ProcLEAR Multifocal has a proven track record of providing outstanding vision, health, and comfort. Now with the addition of an XR range, offering sphere powers up to +/- 20.00D, and ADDs up to +4.00D, even more of your presbyopic patients can enjoy the benefits of ProcLEAR Multifocal.

BALANCED PROGRESSIVE™ TECHNOLOGY— two different, complementary lenses work together to provide exceptional vision at all distances.

PC TECHNOLOGY™ leads to a more comfortable lens-wearing experience—especially important to presbyopic patients whose eyes often grow drier as they age.

SEVEN AVAILABLE ADD POWERS let you more easily fit your presbyopic patients as their condition progresses, offering superior patient retention.




*Based on National Center for Health Statistics data



BALANCED PROGRESSIVE TECHNOLOGY: BETTER VISION BY DESIGN.

Balanced Progressive Technology combines spherical and aspheric optics and unique zone sizes to produce a “D” lens for the dominant eye and an “N” lens for the non-dominant eye. These two lens designs work together to provide a new level of vision quality.



A REVOLUTIONARY MULTIFOCAL SYSTEM.

- Two different, complementary lenses work together to provide exceptional distance, intermediate, and near vision
- Maximizes binocular visual acuity at all ranges of distance with significantly better depth perception and stereopsis

A LENS-FITTING SYSTEM THAT OFFERS MORE CONTROL AND FLEXIBILITY.

- Ability to independently adjust sphere and ADD power by a maximum of +/- 0.50 (ADD can be changed up to +/- 0.50 without affecting sphere power and vice versa)
- A consistent fitting approach lets you keep patients in the same lens design for years to come as their condition progresses

PC TECHNOLOGY: THE SCIENCE OF SUPERIOR COMFORT.

- Lenses remain hydrated—especially important to presbyopic patients whose eyes may be drier due to age
- Only multifocal lens with FDA labeling: “May provide improved comfort for contact lens wearers who experience mild discomfort or symptoms relating to dryness during lens wear.”

PARAMETERS	MATERIAL	WEARING SCHEDULE	WATER CONTENT	BASE CURVE (mm)	DIAMETER (mm)	SPHERE POWER	ADD POWER	LENS DESIGN	HANDLING TINT	LENS CARE	MODALITY
PROCLEAR MULTIFOCAL	omafilcon A	Daily wear	62%	8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	+1.00	D lens N lens	Light blue	Chemical or peroxide disinfection	Monthly replacement
							+1.50				
PROCLEAR MULTIFOCAL XR	omafilcon A	Daily wear	59%	8.7	14.4	+6.25D to +20.00D, -8.50D to -20.00D (0.50D steps after +/- 6.50D)	+2.00	D lens N lens	Light blue	Chemical or peroxide disinfection	Monthly replacement
							+2.50				
						+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+3.00 +3.50 +4.00				
				8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+1.00, +1.50, +2.00, +2.50, +3.00, +3.50, +4.00				

ProcLEAR Multifocal XR is a made-to-order lens and will ship within 10 days from the time of your order.

FOR MULTIFOCAL FITTING CONSULTATION, CALL THE CVSP0 MULTIFOCAL FITTING TEAM AT 800-341-2020 OR VISIT COOPERVISION.COM/MULTIFOCAL.